# USER GUIDE

**Shark**<sup>™</sup>  
**CryoGlow**<sup>™</sup>

The Shark<sup>™</sup> CryoGlow<sup>™</sup> LED mask emits energy in the red and infrared light spectrum for the treatment of fine lines and wrinkles and in the red, blue, and infrared light spectrum for the treatment of mild-to-moderate inflammatory acne.

## TABLE OF CONTENTS

GET TO KNOW YOUR SHARK™ CRYOGLOW™ MASK .....	3
WHAT'S IN THE BOX .....	3
CONTRAINDICATIONS .....	4
WARNINGS & PRECAUTIONS .....	6
LED PARAMETERS .....	8
GUIDANCE & MANUFACTURER'S DECLARATION - EMC .....	8
GETTING STARTED .....	10
FITTING THE MASK .....	11
TIPS FOR BEST RESULTS .....	12
PROGRESS TRACKER .....	13
GET TO KNOW THE REMOTE .....	14
REMOTE SETTINGS .....	15
LED SKIN ROUTINES .....	16
UNDER-EYE COOLING .....	17
LEDS & HOW THEY WORK .....	18
UNDER-EYE CHILL PADS AND HOW THEY WORK .....	19
RUNNING A SESSION .....	20
SELECTING A ROUTINE .....	21
COOLING DURING A ROUTINE .....	22
SELECTING UNDER-EYE COOLING .....	23
SETTINGS, MASK CARE & BATTERY CHARGING .....	24
FAQ .....	25
TROUBLESHOOTING .....	26
ERROR MESSAGES .....	28
ICON EXPLANATIONS .....	29
WARRANTY .....	30
EXPECTED SERVICE LIFE .....	31

SCAN TO GET STARTED Or visit [sharkbeauty.ca](https://sharkbeauty.ca) for more.

# GET TO KNOW YOUR SHARK™ CRYOGLOW™ MASK

### ADJUSTABLE STRAPS

designed for optimal fit and comfort

### MOUTH & NOSE HOLES

for breathability

### SHIELD DESIGN

for even LED coverage

480 LEDS for broad, even coverage

**AUTO-DIM** feature lowers light level when mask is not close to the face

**FOREHEAD PADDING** for a comfortable fit

**SILICONE EYE SHIELDS** for LED light blocking and comfort

**ACTIVE UNDER-EYE COOLING** to soothe under-eyes

## WHAT'S IN THE BOX

Shark™ CryoGlow™ Mask + Straps  
Pre-Assembled

Remote Control with LCD Screen + Control Buttons

Quick Start Guide + User Guide

Universal USB-C Charging Cable + Charging Brick

Clip-on Chill Pads

Storage Bag\*

## CONTRAINDICATIONS

### Do not use the Shark™ CryoGlow™ Mask if you are taking:

- Medication or supplements that cause sensitivity to light (photosensitivity), for example, medications for severe inflammatory acne (i.e. isotretinoin or doxycycline).

### Immediately stop using the device if you experience any discomfort.

### Do not use the Shark™ CryoGlow Mask if you have any of the following medical conditions:

- Retinal and/or macular disease, ocular albinism, or any eye condition that causes sensitivity to light
- Pregnant, nursing, or maybe pregnant
- Pacemaker or other electronic implanted device
- Implants in the face
- Lack of sensation (e.g., numbness)
- Light sensitive skin disorder (e.g., lupus, porphyria)
- Extreme light sensitivity of the eyes (i.e., photophobia)
- Light-induced headaches
- Albinism
- Melasma (i.e., condition that causes patches and spots)
- Sensitivity or allergy to plastic, aluminum, or self-fastening material

### Do not use the Shark™ CryoGlow Mask if you have any of the following conditions on the face, hairline, ears, or neck:

- Recent facial surgery (until you have been cleared by your doctor)
- Suspicious lesions or skin cancer
- Eczema
- Psoriasis
- Irregularly shaped, new, or changing brown spots
- Birthmarks
- Broken skin, skin rash, or blisters
- Tissue inflammation
- Skin infection, cold sore, or warts
- Active current outbreak of herpes simplex virus (cold sores)
- History of herpes outbreaks without receiving preventative treatment
- Facial treatment where the skin has not fully healed, including neurotoxin, dermal filler, micro needling, laser, or chemical peel

## CONTRAINDICATIONS, CONTINUED

### Do not apply the Shark™ CryoGlow™ Mask's cooling treatments on the eyeball or upper eyelid.

### Do not use the Shark™ CryoGlow Mask's cooling treatments if you have any of the following medical conditions:

- Cold hypersensitivity
- Cold urticaria (i.e., hives brought on by cold temperatures)
- Circulatory insufficiency

### PRECAUTIONS: Consult a medical professional before using the Shark™ CryoGlow Mask if:

- You have bruising on the face.
- You are experiencing skin sensitivity or irritation with retinol or benzoyl peroxide.

## WARNINGS & PRECAUTIONS

To reduce the risk of fire, thermal and chemical burns, electric shock, serious injury, or death, follow the precautions in this section and the procedures in this manual. The remote control contains a rechargeable lithium-ion battery. If damaged, the battery can overheat, catch fire, or leak a chemical that can burn skin and eyes.

### SUPERVISION AND INTENDED USERS:

- This equipment is not intended for use by persons with reduced physical, sensory, or mental capabilities, unless they have been given supervision or instruction concerning use of the equipment. Please refer to the contraindications before use.
- Do not use if any of the conditions in “Contraindications” (pages 4-5) apply to you.
- Do not leave the equipment unattended when it is powered on.
- This equipment is not intended for use by children. Children should be supervised to ensure they do not play with the equipment.
  - To reduce the risk of entanglement and potential strangulation, keep cords out of reach of children.
  - Small parts can be inhaled or swallowed by children—if any small parts break off or detach from the product, dispose of safely.
- This equipment is only intended to be used by one person. This person is the intended operator of the equipment

### KEEP DRY.

- Keep the remote control dry. If it gets wet, do not charge or use the battery, as it may be damaged. Follow disposal instructions.
- Water and electricity are a dangerous combination!
  - Do not use this equipment in wet surroundings (e.g., in a bathtub, shower, or swimming pool).
  - Do not let water run into the device.
  - It is very important to keep the device dry. Any contact with water could damage the device. Do not use or store the device in the bathroom or near water.
  - Do not charge or use the device if it has been left outdoors.

### CHARGE SAFELY.

- Follow the charger manufacturer’s instructions.
- Only charge the equipment with the provided USB-C cable and charging brick: Dongguan Yinli Electronics, model YLJXA-T050200. Do not use any other USB-C cable or charging brick.
- Check if the voltage indicated on the equipment corresponds to the local mains voltage before you connect the equipment to avoid risk of electrocution or permanent damage to the equipment.
- If the wall socket used to power the equipment has poor connections, the plug will become hot. Make sure you plug the equipment into a properly installed wall socket.
- Do not charge the device in the bathroom or near water.
- Do not charge the remote after contact with water.
- Always position the equipment so the charger is easy to reach while charging.
- Store USB-C cable properly after charging.

### PROTECT THE REMOTE’S BATTERY FROM DAMAGE.

- Do not try to open the remote or replace the battery.
- Do not crush, drop, or damage the remote.
- Stop use or disconnect charging cable if the remote feels unusually hot or emits smoke or a bad smell.
- If liquid leaks from the remote, avoid contact. If liquid touches skin or eyes, wash with large amounts of water and seek medical advice.
- Follow disposal instructions.

## WARNINGS & PRECAUTIONS, CONTINUED

### USE SAFELY.

- The mask limits your vision. Do not wear while doing anything (such as driving) where limited vision poses a risk.
- CAUTION: The light emitted from this device may be harmful to the eyes. Infrared emitted from this device may cause eye irritation. Do not stare at the light source.
- Do not use the device when it has been left outdoors.
- If you are travelling by air, the mask should be carried in carry-on baggage and not checked in. The device must be completely switched off.
- Do not use the mask outdoors. Only use mask indoors under the following environmental conditions:
  - Temperature: between 50°F (10°C) and 86°F (30°C)
  - A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa.
  - An atmospheric pressure range of 700 hPa to 1060 hPa.
- Do not store or charge the mask face down.
- Skin dryness, irritation, or redness may occur when using an LED routine. Apply a facial moisturizer if dryness occurs.
- To avoid the risk of temporary dry skin or hyperpigmentation, do not run back-to-back LED routines.

### CHECK THE PRODUCT FOR DAMAGE BEFORE EACH USE.

- Inspect the inside of the mask and strap for cleanliness, cracks, or sharp edges.
- Check the chill pads for signs of corrosion.
- These conditions can irritate or damage your skin.

### CLEAN PROPERLY.

- Clean the mask with a damp cloth.
- Avoid using chemical or alcohol-based cleaners.

### DO NOT REPAIR OR MODIFY.

- Do not attempt to repair or modify the mask or remote.
- Do not crush or dispose of in fire.

### DISPOSAL

- Dispose of the device according to local laws and requirements surrounding lithium-ion battery disposal.

### POTENTIAL ELECTROMAGNETIC INTERFERENCE.

- WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- WARNING: Portable radio frequency communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the mask system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- WARNING: Do not leave mask exposed in direct sunlight.

## TRANSPORT AND STORAGE CONDITIONS

- Only transport and store the device in the following environmental conditions:
  - Keep between -4°F (-20°C) and 113°F (45°C). For storage longer than 3 months, keep between 32°F (0°C) and 77°F (25°C).
  - Keep between 41°F (5°C) to 89.6°F (32°C) at a relative humidity up to 85%, non-condensing.
  - Keep between 95°F (35°C) to 140°F (60°C) at a water vapour pressure up to 50 hPa.
  - Atmospheric pressure range: 700 hPa to 1060 hPa.

## LED PARAMETERS

Mode	Better Aging	Skin Clearing		Skin Sustain
Power (mW/cm <sup>2</sup> )	Red 64 ± 5mW/cm <sup>2</sup> , IR 64 ± 5 mW/cm <sup>2</sup>	<b>Step 1:</b> Blue 64 ± 8mW/cm <sup>2</sup> , IR 64 ± 5 mW/cm <sup>2</sup>	<b>Total mode</b> = 128 mW/cm <sup>2</sup>	Blue 59.8 ± 7mW/cm <sup>2</sup> , Red 17 ± 5mW/cm <sup>2</sup> , IR 51.2 ± 5 mW/cm <sup>2</sup>
	<b>Total mode</b> = 128 mW/cm <sup>2</sup>	<b>Step 2:</b> Blue 55 ± 6mW/cm <sup>2</sup> , Red 73 ± 5 mW/cm <sup>2</sup>	<b>Total mode</b> = 128 mW/cm <sup>2</sup>	<b>Total mode</b> = 128 mW/cm <sup>2</sup>
		<b>Step 3:</b> Red 73 ± 5mW/cm <sup>2</sup> , IR 55 ± 5 mW/cm <sup>2</sup>	<b>Total mode</b> = 128 mW/cm <sup>2</sup>	

## GUIDANCE & MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSION

The device is Risk Group 1 as defined in IEC 62471.

The <b>Shark<sup>™</sup> CryoGlow Mask<sup>™</sup></b> is intended for use in the electromagnetic environment specified below. The customer or the user of The <b>Shark CryoGlow Mask</b> should assure that it is used in such an environment.		
EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions	CISPR 11	Home healthcare environment
RF emissions	CISPR 11	Home healthcare environment
Harmonic emissions	IEC 61000-3-2	
Voltage function / flicker emissions	IEC 61000-3-5	

The <b>Shark CryoGlow Mask</b> is intended for use in the electromagnetic environment specified below. The customer or the user of The <b>Shark CryoGlow Mask</b> should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	Up to 8 kV contact Up to 15 kV air	Up to ± 8 kV contact Up to ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	±2kV for main power lines	±2kV for main power lines	Mains power quality should be that of a typical commercial or home environment.
Surge IEC 61000-4-5	±0.5kV 1.5kV differential mode	±0.5 kV to 1.5kV differential mode	Mains power quality should be that of a typical commercial or home environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle <sup>9)</sup> At 0°, 45°, 90°, 135°, 180°, 225°, 270°, & 315° <sup>9)</sup>  0% UT; 1 cycle and 70% <sup>9)</sup> ; UT; 25/30 cycles <sup>10)</sup> Single phase: at 0°  0% UT; 250/300 cycle <sup>10)</sup>	0% UT; 0.5 cycle <sup>9)</sup> At 0°, 45°, 90°, 135°, 180°, 225°, 270°, & 315° <sup>9)</sup>  0% UT; 1 cycle and 70% <sup>9)</sup> ; UT; 25/30 cycles <sup>10)</sup> Single phase: at 0°  0% UT; 250/300 cycle <sup>10)</sup>	Mains power quality should be that of a typical commercial or home environment. If the user of the <b>Shark CryoGlow Mask</b> requires continued operation during power mains interruptions, it is recommended that the <b>Shark CryoGlow Mask</b> be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz magnetic field) IEC 61000-4-8	30 A/m at 50-60Hz	30 A/m at 50-60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or home environment.

NOTE U<sub>i</sub> is the a.c. mains voltage prior to application of the test level.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF IEC 61000-4-6	150KHz - 80MHz - Test level: 3V & 6V - Modulation: AM 1KHz 80%	- 150KHz - 80MHz - Test level: 3V & 6V - Modulation: AM 1KHz 80%	Portable and mobile RF communications equipment should be used no closer to any part of the <b>Shark CryoGlow Mask</b> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance: $d = \left[ \frac{3.5}{f_1} \right] \sqrt{P}$ $d = \left[ \frac{12}{f_2} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{f_1} \right] \sqrt{P}$ $d = \left[ \frac{7}{f_1} \right] \sqrt{P}$ 80 MHz to 800 MHz    800 MHz to 2.7 GHz
Radiated RF IEC 61000-4-3	80MHz - 2.7GHz - Test level: 10V/m - Modulation: AM 1KHz 80%	- 80MHz - 2.7GHz - Test level: 10V/m - Modulation: AM 1KHz 80%	where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). <sup>9)</sup> Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>9)</sup> . Should be less than the compliance level in each frequency range. <sup>9)</sup> Interference may occur in the vicinity of equipment marked with the following symbol:  NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

- a)** The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- b)** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **Shark CryoGlow Mask** is used exceeds the applicable RF compliance level above, the TheraFace Mask should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **Shark CryoGlow Mask**.
- c)** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS

Recommended separation distances between portable and mobile RF communications equipment and the <b>Shark CryoGlow Mask</b>							
The <b>Shark CryoGlow Mask</b> is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <b>Shark CryoGlow Mask</b> can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <b>Shark CryoGlow Mask</b> as recommended below, according to the maximum output power of the communications equipment							
Rated maximum output of transmitter (W)	Separation distance according to frequency of transmitter (m)						
	150 kHz to 80 MHz outside ISM and amateur radio bands	$d = \left[ \frac{3.5}{f_1} \right] \sqrt{P}$	150 kHz to 80 MHz outside ISM and amateur radio bands	$d = \left[ \frac{12}{f_2} \right] \sqrt{P}$	80 MHz to 800 MHz	$d = \left[ \frac{3.5}{f_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz
0.01	0.12		.020		0.035		0.07
0.1	0.38		0.63		0.11		0.22
1	1.2		2.00		0.35		0.70
10	3.8		6.32		1.10		2.21

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.  
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.  
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## GETTING STARTED

Please read the Warnings, Instructions, and Countraindications before first use.

### Charge fully before first use

Plug small end of USB-C cable into charging port on top of remote.

Plug large end of USB-C into power adapter. Plug power adapter into wall outlet.

Allow remote to charge fully (**this should take about 3 hours**). The remote screen will show your charging progress.

Once charged, unplug power adapter and disconnect USB-C cable from remote.

**NOTE:** Unplug before using face mask.  
Mask will not work when plugged into wall outlet.


### Powering ON/OFF

**ON:** To turn on your mask and remote, **long-press** the remote's dial/center button for 2 seconds.

**OFF:** To turn off, **long-press** the dial/center button.

**AUTOMATIC SHUT-OFF:** If the device is not manually turned off, the mask and remote will automatically shut off after completing a session.

The mask fans will continue to run for 60 seconds while the mask is powering down after a session.

## FITTING THE MASK

**NOTE:** Demo video will run upon first use of the remote—please watch it for further fit and comfort instructions. You can also rewatch this video by heading to the remote's settings.

### Step 1 / Adjust the back strap

Start by adjusting the back strap (with the buckles). To tighten or loosen, **slide the buckles simultaneously** together or apart. This strap may cover the top of your ears—adjust until this feels comfortable.


Tighten

Loosen

### Step 2 / Adjust the top strap

To adjust the mask up or down on your face, loosen or tighten the self-fastening strap on the top of your head. Ensure the strap sits in the center of your head (following a center part).

### Step 3 / Adjust chill pad fit (optional)

The chill pads should sit under your eyes, on top of your cheekbones. If the chill pads have poor or inconsistent contact with your under-eye area, **add the clip-on chill pads to improve fit and comfort.**

Remove the mask and take the additional clip-on chill pads out of the box. They are marked "L" and "R"—make sure the letters are pointing toward the chin of the mask. Align the clip-on chill pads with the original chill pads (that are permanently attached to the mask) and snap them in place. You'll feel a click when they are properly attached.

**NOTE:** LEDs work best when they are not in contact with your skin.  
This is why the mask sits at a slight distance from your face.

### Pro Tips:

- 1 The back strap should sit parallel to the floor, and your eyes should be centered in the eye holes.
- 2 Use one hand to hold the back strap in position and the other hand to adjust the top self-fastening strap until the mask is positioned comfortably and securely on your head.

Please visit  
**sharkbeauty.ca**  
for more helpful  
images and videos.

## TIPS FOR BEST RESULTS

### RESULT TIMING

- Results are achieved through consistent use of the same routine daily for 8 weeks or more.
- After consistent use with the same treatment routine, keep going with that routine to build on your progress. Or switch to another treatment routine to work toward another skin goal.
- If you've achieved the skin goals you were working toward, you can use the Skin Sustain routine to help maintain your skin condition over time.

### WHAT TO EXPECT

- You may experience different skin effects in the process of your skin journey when consistently using an LED treatment.
- For the Better Aging routine, skin dryness, irritation and redness might occur. Apply a rich facial moisturizer after each treatment.
- If you experience any undesirable pigmentation or darkening of the skin, please discontinue use and contact SharkNinja at 800-275-9251.

### TAKE “BEFORE AND AFTER” PICTURES

- Be sure to take “before and after” pictures to document your journey while using the LED skin treatments.

## PROGRESS TRACKER

The remote keeps track of your sessions to help you meet your skin goals.

### YOUR ROUTINE

After you finish a treatment, the total number of sessions completed in that routine will appear on the remote's screen.

A reminder will appear above your prior completed treatment when you return.

### PROGRESS

Select “Your Progress” to see your cumulative progress across routines.

### TRACKING

Total Sessions are the number of sessions you have completed, and Streak is how many sessions you have completed in a row.

## GET TO KNOW THE REMOTE

## USING THE REMOTE

### Changing Brightness

1. Turn **ON** the mask using the remote.
2. Turn the dial to view **Settings** and click the dial button to select.
3. Select **Screen Brightness**, then use the dial to increase or decrease the brightness.
4. Click the dial button to make your selection.
5. After selecting a brightness level, you will be sent back to the main menu.

### Changing Language

1. Turn **ON** the mask using the remote.
2. Turn the dial to view **Settings** and click the dial button to select.
3. Select **Language** and use the dial to select the language of your choice.
4. After selecting a language, you will be sent back to the main menu.

### Reset Progress

1. Turn **ON** the mask using the remote.
2. Turn the dial to view **Settings** and click the dial button to select.
3. Select **Reset Progress**, then click the dial button to select **Confirm**.
4. After confirming, you will be sent back to the main menu.

### Factory Reset

1. Turn **ON** the mask using the remote.
2. Turn the dial to view **Settings** and click the dial button to select.
3. Select **Factory Reset**, then click the dial button to select **Confirm**.
4. After confirming, you will be sent back to the main menu.

**Pro Tip:** You can replay the demo video on the remote and access our support page through settings.

## LED SKIN ROUTINES

The mask offers **3 LED routines** to choose from, whether your goals are treating your skin or maintaining it.

### TREATMENT ROUTINES

The Shark® CryoGlow™ LED mask emits energy in the red and infrared light spectrum for the treatment of fine lines and wrinkles and in the red, blue, and infrared light spectrum for the treatment of mild-to-moderate inflammatory acne. Select *either* the **Better Aging** or **Skin Clearing** routine below and use your selection *exclusively* once daily for a minimum of **8 weeks**. For the treatment of **both fine lines/wrinkles and acne**, we recommend using the **Skin Clearing** mode *exclusively* for 8 weeks. Do not stack modes.

#### Better Aging (6 MINUTES)

Enhances skin's appearance by smoothing wrinkles, tightening sagging areas, and balancing texture for a radiant, even look.

#### Skin Clearing (8 MINUTES)

Helps treat mild to moderate inflammatory acne and soothes skin for a clearer and calmer complexion with an automated 3-phase routine.

### MAINTENANCE ROUTINE

Once you have achieved your target results with a treatment routine, use Skin Sustain **once daily** to help maintain your results.

#### Skin Sustain (4 MINUTES)

Helps to boost brightening and helps maintain complexion with daily use.

**Did you lose your streak? Don't worry!**  
**Readjust your goals and continue with your skincare journey.**

No UV light is emitted from this mask. For more information on LEDs, please see page 18.

## UNDER-EYE COOLING

Use cooling in the morning to help soothe your under-eyes and in the evening to relax.

### LED + COOLING

Under-eye cooling will automatically start with any LED routine and has the option to be turned off. Use the remote's chill button to cycle through 3 temperatures to quickly soothe your under-eye area.

### COOLING ONLY

#### Under-Eye Revive (5, 10 or 15 minutes)

An adjustable under-eye cooling **ONLY** routine (no LEDs) for a soothing and refreshing experience designed for optimal relaxation.

## FIND YOUR PERFECT ROUTINE

Scan the QR code to take our quiz to find the best routine for your skin concerns.

## LEDS & HOW THEY WORK

### LED stands for Light Emitting Diode

Different colour LEDs all emit different wavelengths that are measured in nanometers (nm). So when LEDs interact with skin, each colour penetrates different layers.

### LED SKIN PENETRATION

**NOTE:** Infrared wavelengths are invisible to the human eye.


#### LED WAVELENGTHS & TOPICAL BENEFITS

**Blue**  
415nm

Helps kill acne-causing bacteria to treat breakouts.

**Red & Infrared**  
630-830nm

Helps reduce the appearance of fine lines and wrinkles.

## UNDER-EYE CHILL PADS AND HOW THEY WORK

### The under-eye chill pads use the Peltier Effect.

The Peltier Effect describes the movement of heat-producing electricity through both sides of a surface, making one side feel cool and the other side warm. We employ the combination of the Peltier Effect and fan technology to expel excess heat outward—so the metal portion of the mask feels cool to the touch and your under-eye skin experiences instant cooling.

### INSTACHILL TECHNOLOGY USES THE PELTIER EFFECT.

## RUNNING A SESSION

**Note:** Please read the Warnings, Instructions, and Contraindications before first use. Make sure the mask has been fully charged before first use (this could take up to 3 hours).

### HOW TO FIT SHARK™ CRYOGLOW™ INTO YOUR ROUTINE

Take the time with your first few uses to experiment with how the Shark CryoGlow Mask can best fit into your existing routine. Give it a try in the morning to help you feel fresh for the day. If this doesn't work, try it in the evening to help you wind down and relax when you have more time.

#### Step 1 / Put on the mask

Begin with a clean, bare, and dry face.

#### Step 2 / Adjust the back and top straps (Full instructions on page 11)

Adjust the straps to position mask comfortably on your face. Slide buckles simultaneously to adjust the back strap. To adjust the positioning of the back strap on your head, release the top self-fastening strap and gently pull the back strap so it sits parallel to the floor.

Chill pads should be positioned under the eyes, making light contact with skin. **You may need to attach the optional clip-on chill pads available in the box to achieve the best fit.**

The top strap should be positioned in the center of your head, while the back strap should be parallel to the floor.

#### Step 3 / Power on

Turn on the mask by pressing and holding the remote's dial button for 2 seconds. On first use, a short demo video on the remote screen will show you how to fit the mask and interact with the remote.

#### Step 4 / Select and run a routine

Use the remote to choose a routine and begin by selecting "Start." Under-eye cooling will automatically turn on. Adjust the under-eye chill pads' intensity to your desired level with the remote's chill button.

#### Step 5 / Remove and clean the mask

When your session is complete, remove the mask. Fans will run for 60 seconds to power down the chill pads, then the mask will shut off automatically. Once you have taken off the mask, clean the mask with a damp cloth. Continue with your daily skincare routine.

## SELECTING A ROUTINE

**Step 1 /** Use the dial on the remote to go to Routines and scroll through the different **LED** and **Cooling routines** available.

**Step 2 /** Use the dial to select from one of the 3 available LED routines or 1 cooling-only routine. Select your routine by clicking the dial button. Under-eye cooling will automatically start with an LED routine.

## COOLING DURING A ROUTINE

During a routine, you have options to adjust the under-eye chill level, and to pause or end your routine.

### PAUSE

Click dial button to pause session.

### BACK

Press to go to previous screen.

### CHILL

Press to adjust chill level.

### CHILL LEVELS

Level 3

Level 2

Level 1

Off

Shark

## SELECTING UNDER-EYE COOLING

**Step 1** / Use the dial on the remote to scroll to **Under-Eye Revive**.

### INSTACHILL TREATMENT

Use this cooling function to quickly refresh, and soothe the under-eye area.

**Step 2** / Then use the dial to **scroll through run-time options**. Choose one and click the dial button.

**Step 3** / Use the chill button to select your preferred temperature. Mask will automatically shut off at the end of the selected duration.

### Remote settings

- The remote settings are available to adjust from the main menu.
- **Replay tutorial:** Replay demo video shown at first use.
- **Support:** Scan QR code to access our support page.
- **Brightness:** Adjust the brightness of the remote screen.
- **Language:** Choose your preferred language.
- **Reset Progress:** Option to reset your progress level.
- **Factory Reset:** Option to reset your remote/mask.

### Mask care

- Your mask should be cleaned after each use.
- Make sure the mask is turned off and unplugged from any power source.
- Clean the mask with a damp cloth to wipe the interior and exterior of the mask after each use. The cloth should not be dripping wet, as it could damage the electrical components.  
**Ensure the chill pads are fully dry before adding clip-on pads.**
- Allow mask to dry completely before using again.
- **NOTE:** Do not use chemical or alcohol-based cleansers.

### Battery charging

- The run time of a fully charged battery depends on the types of sessions being run. Select the best LED routine for your skin needs and do not stack routines.
- It should take about 3 hours to fully charge the mask. The mask does not work while charging.
- To make sure your mask is always ready for a full session, we recommend charging it after every 2 or 3 uses.

### What happens if I miss a day?

- Consistency is key for results. However, don't be discouraged—just pick up where you left off.

### Is the mask safe to use?

- **Auto-timer:** Mask will shut off automatically after 60 seconds of inactivity once a session is complete.
- **Temperature regulation:** Cooling temperature of the chill pads is regulated.

### Is the mask safe for the eyes?

- Yes, the mask has features to ensure eye safety with LED use.
- **Auto-Dim:** LED lights automatically dim when the device is not close to the face while on.
- **Silicone eye guards:** Ensures safe use with closed or open eyes.

### Does the mask emit UV light?

- No, UV light is not emitted from in this mask. Blue and red LEDs combined can appear purple but do not mimic the effects of UV light.

### How long should I wear the mask?

- The mask's longest LED session is 8 minutes. The longest cooling-only session is 15 minutes. We recommend using only one LED routine per day. Using multiple LED routines per day could lead to dry skin or hyperpigmentation.

### How often should I use the mask?

- For results, use the same LED routine once a day for a minimum of 8 weeks. To maintain results, continue daily usage with Skin Sustain routine.

### When should I see results?

- Results depend on your skin's starting state. With daily use, you may see improvements in about 8 weeks. See page 18 for more information on the biology of your skin when treated with LEDs.

### How do I travel with the mask?

- Please make sure to use a power adapter usable in the country you are traveling in. The mask uses a lithium-ion battery rated at 2.35Ah, 16.92Wh. If you are travelling by air, the battery should be carried in carry-on baggage and not checked in. The device must be completely switched off.

### Whom can I contact for support?

- Should you have any questions or issues with your mask, please contact us by visiting [support.sharkclean.ca](https://support.sharkclean.ca).

### Is the battery removable?

- No, the battery is not removable or replaceable. Should you have any questions or issues with the battery, please contact us by visiting [support.sharkclean.ca](https://support.sharkclean.ca).

## TROUBLESHOOTING

### THE DEVICE DOES NOT POWER ON

#### Possible cause

#### Corrective action

Battery level low

Charge the device using the provided USB-C cable.

Device still will not power on

Long-press the dial button for 2 seconds. If device still will not turn on, please visit [support.sharkclean.ca](https://support.sharkclean.ca).

### THE BATTERY WILL NOT CHARGE

#### Possible cause

#### Corrective action

Possible defective battery

Please visit [support.sharkclean.ca](https://support.sharkclean.ca).

Defective charging cable

Please visit [support.sharkclean.ca](https://support.sharkclean.ca) or replace with a Shark USB-C cable and adapter.

### THE BATTERY TAKES MORE THAN 3 HOURS TO CHARGE

#### Possible cause

#### Corrective action

Wrong power adapter

Make sure to use a 120V power adapter.

**Report any adverse events or issues to customer support at [support.sharkclean.ca](https://support.sharkclean.ca). We monitor safety and performance continually and address any reported issues to maintain high safety standards.**

## TROUBLESHOOTING

### THE CHILL PADS DO NOT POWER ON

#### Possible cause

#### Corrective action

Battery level low

Battery level too low to power the chill pads. Please charge the device using the provided USB-C cable.

The ambient temperature is too high or too low

Turn off the device and let sit for 90 seconds before turning back on. If issue still occurs, stop use and visit [support.sharkclean.ca](https://support.sharkclean.ca).

Chill pads still do not power on

Please contact Shark Customer Service.

### THE COOLING FANS KEEP RUNNING AFTER THE SESSION ENDS

#### Possible cause

#### Corrective action

Device is powering down

Wait 60 seconds for the device to power down after the end of your cooling session.

### THE COOLING FANS ARE NOT RUNNING

#### Possible cause

#### Corrective action

Cooling has been switched off

If cooling is off or turned off during an LED routine, the fans will automatically stop after 60 seconds. This may be why you don't hear them running during the routine.

The cooling fans should be running if cooling is switched on. You can adjust cooling using the dial or the cooling button. Make sure your remote shows that cooling is on and hasn't been turned off.

### THE CHILL PADS ARE DIGGING INTO MY UNDER-EYE AREA

#### Possible cause

#### Corrective action

The straps are not properly placed/fitted on your head

Revisit the "Fitting the Mask" section of this User Guide for step-by-step instructions. The chill pads should only lightly contact your under-eye area.

## ERROR MESSAGES

Error message
Remote Over Temperature. Let remote cool down before using.
Remote Under Temperature. Let remote warm up before using.
Remote error A3. Call Customer Service.
Charger not recognized.
Remove remote from charger.
Low Battery. Charge before use.
Battery error B1. Call Customer Service.
Remote error B2. Charge the remote before next use.
Charging error C1. Disconnect and reconnect to charge.
Fan error F1. Call Customer Service.
Error A1. Call Customer Service.
Error A2. Call Customer Service.

## ICON EXPLANATIONS

Symbol	Explanation
	Type BF Applied Part
	Class II Equipment
	Direct Current
	Refer to User Guide
<b>IP22</b>	Protection against ingress: IP2X: Protected against solid foreign objects of 12.5 mm Ø and greater IPX2: Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH AAMI ES 60601-1:2005/A2:2021; CAN/CSA-C22.2 No. 60601-1:14(Reaffirmed 2022); IEC 60601-1-11:2015/AMD1:2020; IEC 60601-2-83:2019/AMD1:2022	
	Compliant with Medical Standards

## TWO (2) YEAR LIMITED WARRANTY

The Two (2) Year Limited Warranty applies to purchases of the product ("Original Product") made from authorized retailers of SharkNinja Operating LLC. Warranty coverage applies to the original owner and to the original product only and is not transferable.

SharkNinja warrants that the Original Product and battery shall be free from defects in material and workmanship for a period of two (2) years from the date of purchase when it is used under normal household conditions and maintained according to the requirements outlined in the User's Guide, subject to the following conditions and exclusions:

### What is covered by this warranty?

1. The Original Product and/or accessories deemed defective, in SharkNinja's sole discretion, will be repaired or replaced up to two (2) years from the original purchase date.
2. In the event a replacement unit is issued, the warranty coverage ends six (6) months following the receipt date of the replacement unit or the remainder of the existing warranty, whichever is later. SharkNinja reserves the right to replace the Original Product with one of equal or greater value.

### What is not covered by this warranty?

1. Normal wear and tear of wearable parts (such as accessories) which require regular maintenance and/or replacement to ensure the proper functioning of your Original Product, are not covered by this warranty. Replacement parts are available for purchase at [sharkclean.ca/sharkca/parts\\_accessories](https://support.sharkclean.ca/sharkca/parts_accessories).
2. Any Original Product that has been tampered with or used for commercial purposes.
3. Damage caused by misuse, abuse, negligent handling, failure to perform required maintenance (e.g., cleaning device with chemical or alcohol-based cleaners, dropping mask and/or remote, or immersing device in water), or damage due to mishandling.
4. Consequential and incidental damages.
5. Defects or damages caused by repair persons not authorized by SharkNinja. These defects include damages caused in the process of

- shipping, altering, or repairing the Original Product (or any of its accessories or parts) when the repair is performed by a repair person not authorized by SharkNinja.
6. Damages caused by use of parts or accessories which are not Shark Genuine Components.
7. Original Products purchased, used, or operated outside North America.

### How to get service

If your Original Product fails to operate properly while in use under normal household conditions within the warranty period, visit <https://support.sharkclean.ca> for product care and maintenance self-help. Our Customer Service Specialists are also available at **1-800-798-7398** to assist with product support and warranty service options, including the possibility of upgrading to our VIP warranty service options for select product categories. So we may better assist you, please register your product online at [sharkclean.ca/register/guarantee](https://support.sharkclean.ca/register/guarantee) and have the product on hand when you call. SharkNinja will cover the cost for the customer to send in the unit to us for repair or replacement. A fee of \$20.95 (subject to change) will be charged when SharkNinja ships the repaired or replacement unit.

### How to initiate a warranty claim

You must call **1-800-798-7398** to initiate a warranty claim. Vous aurez besoin du reçu comme preuve d'achat. We also ask that you register your product online at [sharkclean.ca/register/guarantee](https://support.sharkclean.ca/register/guarantee) and have the product on hand when you call, so we may better assist you. A Customer Service Specialist will provide you with return and packing instruction information.

### How state/provincial/territorial law applies

1. This warranty gives you specific legal rights, and you also may have other rights that vary based on your home jurisdiction.
2. Some jurisdictions do not permit the exclusion or limitation of incidental or consequential damages so the above might not apply to you.

**SharkNinja Operating LLC**  
Suite 100, 89 A Street  
Needham, MA 02494

## TWO (2) YEAR LIMITED WARRANTY (Continued)

### FCC WARNINGS

**FCC** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to parts 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

1. Reorient or relocate the receiving antenna.
2. Increase the separation between the equipment and receiver.

3. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
4. Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC rules. Operation is subject to the following conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this device.

### ISED Statement

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

## EXPECTED SERVICE LIFE

Two Years

# ***Shark*** BEAUTY™

**XX/XX/2025**

FW312C\_User\_Guide\_510k\_E\_MP\_Mv7

SHARK BEAUTY and SHARK CRYOGLOW are trademarks of SharkNinja Operating LLC.

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© 2025 SharkNinja Operating LLC.
SHARK BEAUTY et SHARK CRYOGLOW sont des marques de commerce de SharkNinja Operating LLC.
© SharkNinja Operating LLC, 2025.